



Product Datasheet

IκB-ε (Phospho-Ser22) Antibody (orb14913)

Description	IκB-ε (Phospho-Ser22) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	054910
Expiration Date	12 months from date of receipt.



Western blot analysis of IκB-ε(Phospho-Ser22)



Immunohistochemical analysis of formalin-fixed tissue showing IκB-ε (Phospho-Ser22) staining.