



Product Datasheet

VASP (Phospho-Ser157) Antibody (orb14914)

Description	VASP (Phospho-Ser157) Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P50552
Expiration Date	12 months from date of receipt.



Western blot analysis of VASP(Phospho-Ser157)



Immunohistochemical analysis of formalin-fixed tissue