



Product Datasheet

4E-BP1 (Phospho-Thr46) Antibody (orb14921)

Description	4E-BP1 (Phospho-Thr46) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13541
Expiration Date	12 months from date of receipt.



Western blot analysis of 4E-BP1(Phospho-...



Immunohistochemical analysis of formalin...