



Product Datasheet

IRS (Phospho-Ser639) Antibody (orb14928)

| | |
|----------------------------|---|
| Description | IRS (Phospho-Ser639) Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P35568 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical analysis of formalin...



Western blot analysis of IRS-1(Phospho-S...