



Product Datasheet

p44/42 MAP Kinase (Phospho-Tyr204) Antibody (orb14936)


Description	p44/42 MAP Kinase (Phospho-Tyr204) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P27361
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...



Western blot analysis of p44/42 MAP Kina...



Immunohistochemical analysis of formalin...