



Product Datasheet

HSP27 (Phospho-Ser78) Antibody (orb14937)


Description	HSP27 (Phospho-Ser78) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P04792
Expiration Date	12 months from date of receipt.



Western blot analysis of HSP27(Phospho-S...



Immunohistochemical analysis of formalin...



Confocal immunofluorescent analysis of H...