



Product Datasheet

NFκB-p105 (Phospho-Ser932) Antibody (orb14940)

Description	NFκB-p105 (Phospho-Ser932) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P19838
Expiration Date	12 months from date of receipt.



Western blot analysis of NFκB-p105(Phosp...



Immunohistochemical analysis of formalin...