



Product Datasheet

P38 MAPK (Phospho-Thr180) Antibody (orb14941)

Description	P38 MAPK (Phospho-Thr180) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16539
Expiration Date	12 months from date of receipt.



Western blot analysis of P38 MAPK(Phosph...



Immunohistochemical analysis of formalin...