

Product Datasheet

Recombinant Enterokinase, Light Chain, Porcine (orb1494616)

Description	GenScript Enterokinase is a highly purified recombinant porcine enterokinase. The enzyme has been ex...
Preservatives	GenScript enterokinase, formulated using GenScript's proprietary technology, can be shipped at room temperature. It will remain stable at 37°C for one week without losing any activity.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at -20°C after delivery.
Note	For research use only
Application notes	Reconstituted in ddH ₂ O or PBS at 100 µg/ml.
Source	<i>P. pastoris</i>
MW	Theoretical MW: 21,880 Da The apparent MW on SDS-PAGE: about 40,000 Da
Expiration Date	6 months from date of receipt.