

Product Datasheet

Recombinant IP-10/CRG-2/CXCL10, Rat (orb1494635)

Description	C-X-C motif chemokine 10 (CXCL10) also known as interferon gamma-induced protein 10 (IP-10) or small...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Rat IP-10°CRG-2°CXCL10 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Rat IP-10°CRG-2°CXCL10 should be stable up to 1 week at 4°C or up to 3 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	IPLARTVRCTCIDFHEQPLRPRAIGKLEIIPASLSCPHVEIIATMKKNNEKRCLNPES EA IKSLKAVSQRRSKRAP
Source	HEK 293
MW	8.7 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.