

Product Datasheet

Recombinant IL-1 RA, Human(HEK 293-expressed) (orb1494643)

Description	IL-1 Receptor Antagonist, also known as IL-1RA, ICIL-1RA, IRAP and IL-1RN, is a member of the interl...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant human IL-1 Receptor Antagonist remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Human IL-1 Receptor Antagonist should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	RPSGRKSSKMQAFRIWDVNQKTFYLRNNQLVAGYLQGPNNLEEKIDVVPIEPH ALFLGIHGGKMCLSCVKSGDETRLQL EAVNITDLSNRKQDKRFAFIRSDSGPTTSFESAACPGWFLCTAMEADQPVSLTN MPDEGVMVTKFYFQEDE
Source	HEK 293
MW	18-23 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.