

Product Datasheet

Recombinant GRO- α /KC/CINC-1/CXCL1, Rat (orb1494647)

Description	GRO/MGSA/CXCL1 has chemotactic activity for neutrophils. It may play a role in inflammation and exer...
Endotoxins	< 0.2 EU/ μ g, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Rat GRO/MGSA/CXCL1 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Rat GRO/MGSA/CXCL1 should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH ₂ O or PBS at 100 μ g/ml.
Protein Sequence	APVANELRCQCLQTVAGIHFKNIQSLKVMPPGPHCTQTEVIATLKNGREACLDPE APMVQKIVQKMLKGVPK
Source	HEK 293
MW	~7.8 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.