

Product Datasheet

Recombinant VEGF-D, Human (orb1494659)

Description	Vascular Endothelial Growth Factor (VEGF)-D, also known as c-Fos-induced growth factor (FIGF), is a ...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Human VEGF-D remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Human VEGF-D should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH ₂ O or PBS at 100 μg/ml.
Protein Sequence	FAATFYDIETLKVIDEWQRTQCSPRETCVEVASELGKSTNTFFKPPCVNVFRCG GCCNEESLICMNTSTSYISKQLFEI SVPLTSVPELVPVKVANHTGCKCLPTAPRHPYSIIRR
Source	CHO
MW	18-19 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.