

Product Datasheet

Recombinant TNF R I, Human (orb1494664)

Description	TNF Receptor Type I, is also known as TNF R-p55/p60 and TNFRSF1A. It is a type I transmembrane prote...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Human TNF Receptor Type I remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Human TNF Receptor Type I should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	DSVCPQGKYIHPQNN SICCTKCHKGTYLYNDCPGPGQDTCRECESGSFTASEN HLRHCLSCSKCRKEMGQVEISSCTVD RDTVCGCRKNQYRHYWSENLFQCFNCSLCLNGTVHLSCQEKQNTVCTCHAGFF LRENECVSCSNCKKSLECTKLCLPQIE N
Source	CHO
MW	28~35 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.