

Product Datasheet

Recombinant M-CSF, Rat (orb1494682)

Description	Macrophage-Colony Stimulating Factor (M-CSF), also known as Colony Stimulating Factor-1 (CSF-1), can...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Rat Macrophage-Colony Stimulating Factor (M-CSF) remains stable up to 6 months at -80 °C from date of receipt. Upon reconstitution, rrM-CSF should be stable up to 1 week at 4 °C or up to 2 months at -20 °C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	EVSEHCSTMIGNHGLQILQLIDSQMETACLIIEYKFVDQEQLDDPVCYLKKAFLVQVIEETMRFKDNTPNANATERLQ ELSMKLNSCFIKDYKEQNEACVQTYKESPLRLLEKIKNFFNETKNFLEKDOWNIFSK NCNDSLAKCSSRDVVTKP
Source	CHO
MW	32-40 kDa, observed by non-reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.