

Product Datasheet

Recombinant MCP-3/CCL7, Human(CHO-expressed) (orb1494687)

Description	Chemokine (C-C motif) ligand 7 (CCL7) is a small cytokine that was previously called monocyte-specif...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant human MCP-3°CCL7 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, human MCP-3°CCL7 should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	QPVGINTSTTCCYRFINKKIPKQRLESYRRRTTSSHCPREAVIFKTKLDKEICADPTQ KVVQDFMKHLDKKTQTPKL
Source	CHO
MW	10-18 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.