

Product Datasheet

Recombinant IL-6, Human(CHO-expressed) (orb1494697)

Description	Interleukin-6 (IL-6), also known as BSF-2, CDF and IFNB2, is a pleiotropic cytokine that acts in bot...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Human Interleukin-6 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Human Interleukin-6 should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	PVPPGEDSKDVAAPHRQPLTSSERIDKQIRYILDGISALRKETCNKSNMCESSKEA LAENNLNLPKMAEKDGCFFQSGFNE ETCLVKIITGLLEFVYLEYLQNRFSSEEQARAVQMSTKVLIQFLQKAKNLDAIT TPDPTTNASLLTKLQAQNQLQD MTTHLILRSFKEFLQSSLRALRQM
Source	CHO
MW	21-23 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.