

Product Datasheet

Recombinant IL-3, His, Rat (orb1494708)

Description	Interleukin-3 (IL-3), also known as MCGF, Multi-CSF, HCGF and P-cell stimulation factor, belongs tot...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant ratInterleukin-3 (IL-3), His remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rat Interleukin-3 (IL-3), His should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	ISDRGSDAHLLRLDCRTIALEILVKLPYPQVSGLNNSDDKANLRNSTLRRVNL DEFLKSQEEFDSQDTTDDIKSKLQKL KCCIPAAASDSVLPGVYNKDLDDFKKKLRFYVIHLKDLQPVSVSRPPQPTSSSDN FRPMTVECHHHHHH
Source	CHO
MW	22-34 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.