

Product Datasheet

Recombinant IL-17A, His, Human (orb1494715)

Description	Interleukin-17A (IL-17A), also known as CTLA-8 and IL-17, is a proinflammatory cytokine belonging to ...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant human Interleukin-17A (IL-17A), His remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, human Interleukin-17A (IL-17A), His should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH ₂ O or PBS at 100 μg/ml.
Protein Sequence	GITIPRNP GCPNSEDKNFPRTVMVNLNIHNRNTNTNPKRSSDYNRSTSPWNLH RNEDPERYP SVIWEAKCRHLGCINAD GNVDYHMNSVPIQQEILVLRREPPHCPNSFRLEKILVSVGCTCVTPIVHHVAHHH HHH
Source	CHO
MW	14-22 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.