

Product Datasheet

Recombinant IL-10, Human(CHO-expressed) (orb1494723)

Description	Interleukin-10 (IL-10), initially known as Cytokine Synthesis Inhibitory Factor (CSIF), belongs to t...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Human Interleukin-10 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Human Interleukin-10 should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	SPQGQTQSENSCTHFPGNLPNMLRDLRDAFSRVKTFQMKDQLDNLKESLLE DFKGYLGCQALSEMIQFYLEEVMPQA ENQDPDIKAHVNSLGENLKTLLRRLRRCHRFLPCENKSKAVEQVKNAFNKLQEK GIYKAMSEFDIFINYIEAYMTMKIRN
Source	CHO
MW	18-19 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.