

Product Datasheet

Recombinant IFN- γ , Rat (CHO-expressed) (orb1494724)

Description	Interferon- γ (IFN- γ), also known as Type II interferon or immune interferon, is a cytokine produced...
Endotoxins	< 0.2 EU/ μ g, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Rat Interferon gamma (IFN- γ) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rIFN- γ should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH ₂ O or PBS at 100 μ g/ml.
Protein Sequence	QGTLIESLESKKNYFNSSSSMDAMEGKSLLLDIWRNWQKDGNTKILESQIISFYLRLL FEVLKDNQAISNNISVIESHLITN FFSNSKAKKDAFMSIAKFEVNNPQIQHKAVNELIRVIHQLSPESSLRKRKRSRC
Source	CHO
MW	15-25 kDa, observed by non-reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.