

Product Datasheet

Recombinant GRO- γ /CXCL3, Human(CHO-expressed) (orb1494731)

Description	Chemokine (C-X-C motif) ligand 3 (CXCL3) is a small cytokine belonging to the CXC chemokine family t...
Endotoxins	< 0.2 EU/ μ g, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Human GRO-gamma/CXCL3 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Human GRO-gamma/CXCL3 should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μ g/ml.
Protein Sequence	ASVVELRCQCLQTLQGIHLKNIQSVNVRSPGPHCAQTEVIATLKNKGKACLNPA SPMVQ KIIEKILNKGSTN
Source	CHO
MW	9 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.