

## Product Datasheet

### Recombinant Flt-3L, Human (orb1494745)

<b>Description</b>	Fms-related tyrosine kinase 3 ligand (Flt3L) is growth factor stimulates the proliferation and differ...
<b>Endotoxins</b>	< 0.2 EU/μg, determined by LAL method.
<b>Preservatives</b>	Lyophilized after extensive dialysis against PBS.
<b>Form/Appearance</b>	Sterile Filtered White lyophilized (freeze-dried) powder.
<b>Storage</b>	Lyophilized recombinant human Fms-related tyrosine kinase 3 ligand (Flt-3 Ligand) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rhFlt-3L should be stable up to 1 week at 4°C or up to 2 months at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Reconstituted in ddH2O or PBS at 100 μg/ml.
<b>Protein Sequence</b>	TQDCSFQHSPISDFAVKIRELSDYLLQDYPVTVASNLQDEELCGGLWRLVLAQR WMERLKTVAGSKMQGLLERVNTIEH FVTKCAFQPPPSCLRFVQTNISRLQETSEQLVALKPWITRQNFRCLELQCQPD SSTLPPPWSRPLEATAPTA
<b>Source</b>	CHO
<b>MW</b>	Multiple bands between 18~22 kDa, observed by reducing SDS-PAGE.
<b>Expiration Date</b>	6 months from date of receipt.