

Product Datasheet

Recombinant Ck beta 8-1/CCL23, Human (orb1494757)

Description	Chemokine (C-C motif) ligand 23 (CCL23) is a small cytokine belonging to the CC chemokine family tha...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant human Ck beta 8-1° CCL23 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, human Ck beta 8-1° CCL23 should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Protein Sequence	MLWRRKIGPQMTLSHAAGFHATSADCCISYTPRSIPCSLLESYFETNSECSKPGVI FLTKKGRRFCANPSDKQVQVCVRM LKLDTRIKTRKN
Source	CHO
MW	10.5 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.