

## Product Datasheet

### Recombinant Streptavidin (orb1494771)

<b>Description</b>	Streptavidin is a biotin-binding protein found in the culture broth of the bacterium <i>Streptomyces av...</i>
<b>Form/Appearance</b>	Sterile Filtered White lyophilized (freeze-dried) powder.
<b>Storage</b>	The powder may be stored at 4°C for several weeks. For long-term storage, store at -20°C.
<b>Note</b>	For research use only
<b>Source</b>	<i>Escherichia coli</i> .
<b>MW</b>	13,500 Da
<b>Expiration Date</b>	6 months from date of receipt.