

Product Datasheet

Recombinant IL-15, Human (orb1494832)

Description	Interleukin-15 (IL-15) is a cytokine with structural similarity to IL-2. Like IL-2, IL-15 binds to a...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant Human Interleukin-15 (rhIL-15) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rhIL-15 should be stable up to 1 week at 4°C or up to 3 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH ₂ O or PBS at 100 μg/ml.
Protein Sequence	MNWWNVISDL KKIEDLIQSM HIDATLYTES DVHPSCKVTA MKCFLELQV ISLESGDASI HDTVENLIL ANNSLSSNGN VTESGCKECE ELEEKNIKEF LQSFVHIVQM FINTS
Source	Escherichia coli.
MW	12.8 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.