

Product Datasheet

Recombinant IL-15, His, Human (orb1494833)

| | |
|--------------------------|--|
| Description | Interleukin-15 (IL-15) is a cytokine with structural similarity to IL-2. Like IL-2, IL-15 binds to a... |
| Endotoxins | < 0.2 EU/μg, determined by LAL method. |
| Preservatives | Lyophilized after extensive dialysis against PBS. |
| Form/Appearance | Sterile Filtered White lyophilized (freeze-dried) powder. |
| Storage | Lyophilized recombinant Human IL-15(His-tag) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Human IL-15(His-tag) should be stable up to 1 week at 4°C or up to 3 months at -20°C. |
| Note | For research use only |
| Application notes | Reconstituted in ddH ₂ O or PBS at 100 μg/ml. |
| Protein Sequence | MHHHHHHNWW NVISDLKKIE DLIQSMHIDA TLYTESDVHP SCKVTAMKCF LLELQVISLE SGDASIHDTV ENLILANNS LSSNGNVTES GCKECEELEE KNIKEFLQSF VHIVQMFINT S |
| Source | Escherichia coli. |
| MW | 13.6 kDa, observed by reducing SDS-PAGE. |
| Expiration Date | 6 months from date of receipt. |