

Product Datasheet

Recombinant IL-10, Rat (orb1494834)

Description	Interleukin-10 (IL-10), also known as cytokine synthesis inhibitory factor (CSIF), is an anti-inflam...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant rat Interleukin-10 (IL-10) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, IL-10 should be stable up to 2 weeks at 4°C or up to 3 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O at 100 μg/ml.
Protein Sequence	MSKGHSIKG DNNCTHFPVS QTHMLRELRA AFSQVKTFQ KKDQLDNIVL TDSLLQDFKG YLGCQALSEM IKFYLVEMP QAENHGPEIK EHLNSLGEKL KTLWIQLRRC HRFLPCENKS KAVEQVKNDK NKLQDKGVYK AMNEFDIFIN CIEAYVTLKM KN
Source	Escherichia coli.
MW	18.7 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.