

Product Datasheet

Recombinant IFN- γ , Rat (orb1494841)

| | |
|--------------------------|--|
| Description | Interferon gamma (IFN- γ), also known as Type II interferon, is a cytokine produced primarily by T-ly... |
| Endotoxins | < 0.2 EU/ μ g, determined by LAL method. |
| Preservatives | Lyophilized after extensive dialysis against PBS. |
| Form/Appearance | Sterile Filtered White lyophilized (freeze-dried) powder. |
| Storage | Lyophilized recombinant Rat IFN- γ remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, recombinant Rat IFN- γ should be stable up to 1 week at 4°C or up to 2 months at -20°C. |
| Note | For research use only |
| Application notes | Reconstituted in ddH ₂ O or PBS at 100 μ g/ml. |
| Protein Sequence | GT LIESLES LKN YFNSSSMDAM EGKSLLLDIW RNWQKDGNTK ILESQIISFY LRLFEVLKDN QAISNNISVI ESHLITNFFS NSKAKKDAFM SIAKFEVNNP QIQHKAVNEL IRVIHQLSPE SSLRKRKRSR C |
| Source | Escherichia coli. |
| MW | 15.5 kDa, observed by reducing SDS-PAGE. |
| Expiration Date | 6 months from date of receipt. |