

## Product Datasheet

# Recombinant Glutathione S-Transferase (orb1494857)

<b>Description</b>	Glutathione S-Transferase (GST), an antioxidant enzyme, is involved in the primary cellular defense ...
<b>Form/Appearance</b>	Sterile Filtered White lyophilized (freeze-dried) powder.
<b>Storage</b>	This product is lyophilized in 1X PBS, and it remains stable for several weeks if stored at 4°C. It may be stored it at -20°C for longer periods.
<b>Note</b>	For research use only
<b>Application notes</b>	Add 1 ml of dH2O to make the final concentration of 1 mg/ml in 1X PBS
<b>Protein Sequence</b>	MSPILGYWKI KGLVQPTRLL LEYLEEKYEE HLYERDEGDK WRNKKFELGL EFPNLPYYID GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV DFLSKLPEML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK KRIEAIPOID KYLKSSKYIA WPLQGWQATF GGGDHPPKSD LVPRGSPGIH RD
<b>Source</b>	Escherichia coli.
<b>MW</b>	26.98kD
<b>Expiration Date</b>	6 months from date of receipt.