

## Product Datasheet

### Recombinant FGF-9, Human (orb1494872)

<b>Description</b>	Fibroblast Growth Factor-9 (FGF-9) is a heparin binding growth factor that belongs to the fibroblast...
<b>Endotoxins</b>	< 0.2 EU/μg, determined by LAL method.
<b>Preservatives</b>	Lyophilized after extensive dialysis against PBS.
<b>Form/Appearance</b>	Sterile Filtered White lyophilized (freeze-dried) powder.
<b>Storage</b>	Lyophilized recombinant Human Fibroblast Growth Factor-9(rhFGF-9) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rhFGF-9 should be stable up to 2 weeks at 4°C or up to 3 months at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Reconstituted in ddH2O at 100 μg/ml.
<b>Protein Sequence</b>	MAPLGEVGN YFGVQDAVPF GNVPVLPVDS PVLLSDHLGQ SEAGGLPRGP AVTDL DHLKG ILRRRQLYCR TGFHLEIFPN GTIQGTRKDH SRFGILEFIS IAVGLVSIRG VDSGLYLGMM EKGELYGSEK LTQECVFREQ FEENWYNTYS SNLYKHVDTG RRYVALNKD GTPREGTRTK RHQKFTHFLP RPVDPDKVPE LYKDILSQS
<b>Source</b>	Escherichia coli.
<b>MW</b>	23.4 kDa, observed by reducing SDS-PAGE.
<b>Expiration Date</b>	6 months from date of receipt.