

## Product Datasheet

### Recombinant FGF-8e, Human (orb1494874)

<b>Description</b>	Fibroblast Growth Factor 8e (FGF-8e) is a cytokine belonging to the heparin-binding FGF family, whic...
<b>Endotoxins</b>	< 0.2 EU/μg, determined by LAL method.
<b>Preservatives</b>	Lyophilized after extensive dialysis against PBS.
<b>Form/Appearance</b>	Sterile Filtered White lyophilized (freeze-dried) powder.
<b>Storage</b>	Lyophilized recombinant Human FGF 8e remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, Human FGF 8e should be stable up to 1 week at 4°C or up to 3 months at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Reconstituted in ddH <sub>2</sub> O or PBS at 100 μg/ml.
<b>Protein Sequence</b>	MQEGPGRGPA LGRELASLFR AGREPQGVSQ QHVREQLSVT DQLSRRLIRT YQLYSRTSGK HVQVLANKRI NAMAEDGDPF AKLIVETDTF GSRVRVRGAE TGLYICMNKK GKLIKSNKG GKDCVFTEIV LENNYTALQN AKYEGWYMAF TRKGRPRKGS KTRQHOREVH FMKRLPRGHH TTEQSLRFEF LNYPFTRSL RGSQRTWAPE PR
<b>Source</b>	Escherichia coli.
<b>MW</b>	24.3 kDa, observed by reducing SDS-PAGE.
<b>Expiration Date</b>	6 months from date of receipt.