

Product Datasheet

Recombinant Epigen, Human (orb1494889)

Description	Epigen is a cytokine belonging to the Epidermal Growth Factor (EGF) superfamily, which also includes...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant human Epigen (rhEpigen) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rhEpigen remains stable up to 2 weeks at 4°C or up to 3 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH2O at 100 μg/mL.
Protein Sequence	MAVTVTPPI TAQQADNIEG PIALKFSHLC LEDHNSYCIN GACAFHHELE KAICRCFTGY TGERCEHLTL TSYA
Source	Escherichia coli.
MW	8.1 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.