

Product Datasheet

Recombinant DKK-1, Human (orb1494896)

Description	Dickkopf related protein 1 (DKK1) is a chemokine that belongs to the DKK protein family, which also...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Lyophilized after extensive dialysis against PBS.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Lyophilized recombinant human DKK1 (rhDKK1) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rhDKK1 should be stable up to 1 week at 4°C or up to 2 months at -20°C.
Note	For research use only
Application notes	Reconstituted in ddH ₂ O or PBS at 100 μg/ml.
Protein Sequence	TLNSVLNSNAIKNLPPPLGGAAGHPGSAVSAAPGILYPGGNKYQTIDNYQPYP CAEDEC GTDEYCASPTRGGDAGVQICLACRKRRCMRHAMCCPGNYCKNGICVSSDQ NHFRGEIE ETITESFGNDHSTLDGYSRRITLSSKMYHTKGQEGSVCLRSSDCASGLCCARHF WSKICK PVLKEGQVCTKHRRKGGSHGLEIFQRCYCGEGLSRIQKDHHQASNSSRLHTCQR H
Source	Escherichia coli.
MW	17-22 kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.