

Product Datasheet

Recombinant 3C protease, His, Human (orb1494911)

Description	Recombinant PreScission Protease (rPSP) from human rhinovirus (HRV 3C) is a recombinant restriction-...
Endotoxins	< 0.2 EU/μg, determined by LAL method.
Preservatives	Sterile liquid solution contains 50 mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 50% Glycerol, pH 7.5.
Form/Appearance	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Recombinant PreScission Protease (rPSP) remains stable up to 1 year at -20°C from date of receipt. Please avoid freeze-thaw cycles.
Note	For research use only
Application notes	Reconstituted in ddH2O or PBS at 100 μg/ml.
Source	Escherichia coli.
MW	24kDa, observed by reducing SDS-PAGE.
Expiration Date	6 months from date of receipt.