

---

## Product Datasheet

### Murine Interleukin-6 (rMuIL-6) (orb1495008)

**Description**

Interleukin 6 is a multi functional protein

defense, acute phase ...

**Endotoxins**

Less than 1EU/mg of rMUL-6 as determined by LAL method.

**Preservatives**

Lyophilized from a 0.2mm filtered concentrated solution in PBS, pH 7.4.

**Form/Appearance**

Sterile Filtered White lyophilized (freeze-dried) powder.

**Storage**

This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

**Note**

For research use only

**Application notes**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml. Stock solutions should be apportioned into working aliquots and stored at -20°C. Further dilutions should be made in appropriate buffered solutions.

**Protein Sequence**

FPTSQVRRGD FTEDTTPNRP VYTTSQVGGL  
ITHVLWEIVE MRKELCNGNS DCMNDDALA  
ENNLKLPFIQ RNDGCYQTGY NQEICLLKIS  
SGLLEYHSYL EYMKNLKDKN KKDKARVLQR  
DTETLIHIFN QEVKDLHKIV LPTPISNALL  
TDKLESQKEW LRTKTIQFIL KSLEEFKVT  
LRSTRQT

**Source**

Escherichia coli.

**MW**

Approximately 21.7 kDa, a single non-glycosylated polypeptide chain containing 187 amino acids

**Expiration Date**

6 months from date of receipt.