
Product Datasheet

Mouse IL-1 Alpha (rMuIL-1 alpha) (orb1495016)

Description

IL-1 alpha is a non-secreted

variety of cells including monoc...

Endotoxins

Less than 1EU/mg of rMUL-1 alpha as determined by LAL method.

Preservatives

Lyophilized from a 0.2mm filtered concentrated solution in PBS, pH 7.4.

Form/Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Storage

This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Note

For research use only

Application notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at -20°C. Further dilutions should be made in appropriate buffered solutions.

Protein Sequence

SAPYTYQSDL RYKLMKLVQR
 KFVMNDSLQ TIYQDVDKHY
 LSTTWLNDLQ QEVKFDMYAY
 SSGGDDSKYP VTLKISDSQL FVSAQGEDQP
 VLLKELPETP KLITGETDL IFFWKSINSK
 NYFTSAAYPE LFIATKEQSR VHLARGLPSM
 TDFQIS

Source

Escherichia col

MW

Approximately 17.9 kDa, a single non-glycosylated polypeptide chain containing 156 amino acids.

Expiration Date

6 months from date of receipt.