
Product Datasheet

**Human Platelet-derived growth factor -BB
(rHuPDGF-BB) (orb1495024)**

Description

Human Platelet-derived growth

factor (PDGF) presenting in serum b...

Endotoxins

Less than 1EU/μg of rHuPDGF-BB as determined by LAL method.

Preservatives

Lyophilized from a 0.2μm filtered concentrated solution in PBS, pH 7.4.

Form/Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Storage

This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Note

For research use only

Application notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at -20°C. Further dilutions should be made in appropriate buffered solutions.

Protein Sequence

SLGSLTIAEP AMIAECKTRT
 EVFEISRRLI DRTNANFLVW
 PPCVEVQRCS GCCNNRNVQC
 RPTQVQLRPV QVRKIEIVRK
 KPIFKKATVT LEDHLACKCE
 TVAAARPVT

Source

Escherichia coli.

MW

Approximately 12.4 KDa, a disulfide-linked homodimeric protein containing two 109 amino acid residues polypeptide (B