

Product Datasheet

β -Actin Recombinant Protein (orb1495086)

Description	β -Actin Recombinant Protein
Species/Host	E. coli
Preservatives	PBS, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only
Protein Sequence	MDDDIAALVVDNGSGMCKAGFAGDDAPRAVFPSIVGRPRRHQGMVGMGQKD SYVGDEAQSKRGILTLKYPIEHGIVTNWDDMEKIWHHTFYNELRVAPEEHPVLLT EAPLNPKANREKMTQIMFETFNTPAMYVAIQAVLSLYASGRRTTGIVMDSGDGVT HTVPIYEGYALPHAILRLDLAGRDLTDYLMKILTERGYSFTTTAEREIVRDIKEKLCY VALDFEQEMATAASSSSLEKSYELPDGQVITIGNERFRCPEALFQPSFLGMESCGI HETTFNSIMKCDVDIRKDLYANTVLSGGTTMYPGIADRMQKEITALAPSTMKIKIIA PPERKYSVWIGGSILASLSTFQQMWISKQEYDESGPSIVHRKCF
MW	~43kDa
Uniprot ID	P60709
Solubility (25°C)	PBS, PH7.4
Expiration Date	6 months from date of receipt.