

Product Datasheet

NDRG1 Recombinant Protein (orb1495088)

| | |
|--------------------------|--|
| Description | NDRG1 Recombinant Protein |
| Species/Host | E. coli |
| Preservatives | PBS, pH 7.4 |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Tag | His-tag |
| Note | For research use only |
| Protein Sequence | HMSREMQDVDLAEVKPLVEKGETITGLLQEFDVQEQDIETLHGSVHVTLCGTPK GNRPVILTYHDIGMNHKTCYNPLFNIEDMQEITQHFAVCHVDAPGQDGAASFP AGYMYPSMDQLAEMLPGLVQQFGLKSIIGMGTGAGAYILTRFALNNPEMVEGLV LINVNPCAEGWMDWAASKISGWTQALPDMVVSHLFGKEEMQSNVEVVHTYRQ HIVNDMNPGLHLFINAYNSRRDLEIERPMPGTHVTTLQCPALLVVGDS SPAVDA VVECNSKLDPTKTTLLKMADCGGLPQISQPAKLAEAFKYFVQGMGYMPSASMTR LMRSRTASGSSVTSLDGTRRSRSHSTSEGTRRSRSHSTSEGTRRSRSHSTSEGAHLDTIP NSGAAGNSAGPKSMEVSCLE |
| MW | ~43kDa |
| Uniprot ID | Q92597 |
| Solubility (25°C) | PBS, PH7.4 |
| Expiration Date | 6 months from date of receipt. |