

Product Datasheet

A44R Recombinant Protein (orb1495091)

Description	A44R Recombinant Protein
Species/Host	E. coli
Preservatives	2M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only
Protein Sequence	HMDKIKITIDSKIGNVVTISYNLEKITIDVTPKKKKEKDVLLAQSVAVEEAKDVKVEE KNIIDIEDDDMDIENTLE
MW	~18kDa
Uniprot ID	Q3I957
Solubility (25°C)	2M Urea, PH7.4
Expiration Date	6 months from date of receipt.