

Product Datasheet

A27L Recombinant Protein (orb1495094)

Description	A27L Recombinant Protein
Species/Host	E. coli
Preservatives	0.5M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only
Protein Sequence	HMNKLFSVFSTRYDAKAYLYTEYFLFLNINHYLLSGSDLFINVAYGAVSFSSPISVP DYIMEALTFKACDHIMKSGDLKYTYAFTKKVKDLFNTKSDSVYQYVRLHEMSYDG VSEDTDDDDEVFAILNLSIDSSVDYRNRVLLLTPEVASLRKEYSETEPDYKYLMD EEVPAYDKHLSKPITNTGIEEPHATGGDKEDQPIKVVHPPPNNDKDDAIKYPNPLE DPNYVPTNTRTAIGIADYQLVINKLIEWLDKCEEECGNDGEFKTDLEEAKRKLTEL NEELSDKLSKIRTLERDSVYKTERIDRLTTEIKELRDMQQNGTDDGSDSSEIDKKT IRELRESLDMEREMRSELEKELDTIRDGKVDGSCQRELELSRMWLKQRDDDLRA EIDKRRNVEWELSRLRMDIKECDKYKEDLDKDKTTISTYVSRISTLESEIAKYQQD RDTLVSVHRELEERRHVRDLESRLDECTRNQEDTQEVDALRSRIRELENKLT CIESGGGNLTEISRLQSRISDLERQLSECRGNATEITILE
MW	~64kDa
Uniprot ID	O02685
Solubility (25°C)	0.5M Urea, PH7.4
Expiration Date	6 months from date of receipt.