

Product Datasheet

B2R Recombinant Protein (orb1495096)

Description	B2R Recombinant Protein
Species/Host	E. coli
Preservatives	PBS, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag;Sumo-Tag
Note	For research use only
Protein Sequence	HMTQLPILLLLISLVYATPSPQTSKKIGDDATISCSRNTNYYVMSAWYKEPNSII LLAAKSDVLYFDNYTKDKISYDSPYDDLVTITIKSLTAGDAGTYICAFFMTSTTND TDKVDYEEYSIELIVNTDSESTIDIILSGSTPETISEKPEDIDNSNCSSVFEIATPEPIT DNVEDHTDVTYTSDSINTVNASSGESTTDETPPEITDKEDHTVTDVSYTTVS TSSGIVTTKSTDDADLYDYNNDTVPTTVGGSTTSISNYKTKDFVEIFGITLLII LSAVAIFCITYYICNKHPRKYKTENKVLE
MW	~48kDa
Uniprot ID	Q8AZ47
Solubility (25°C)	PBS, PH7.4
Expiration Date	6 months from date of receipt.