

Product Datasheet

C19L Recombinant Protein (orb1495097)

Description	C19L Recombinant Protein
Species/Host	E. coli
Preservatives	2M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag;Sumo-Tag
Note	For research use only
Protein Sequence	HMWPFASVPAGAKCRLVETLPENMDFRSDHLTTFECFNEIITLAKKYIYIASFCCN PLSTTRGALIFDKLKEVSEKGIKIIVLLDERGKRNLGELQSHSPDINFITVNIKKNN VGLLLGCFWVSDDERCYVGNASFTGGSIHITIKTLGVYSDYPPLATDLRRRFDTFK AFNSAKNSWLNLCSAACCLPVSTAYHIKNPIGGVFFTDSPHELLGYSRDLDTDVV IDKLKSAKTSIDIEHLAIVPTTRVDGNSYYWPDYNSIIEAAINRGVKIRLLVGNWDK NDVYSMATARSLDALCVQNDLSVKVFTIQNNTKLLIVDDEYVHITSANFDGTHYQ NHGFVSFNSIDKQLVSEAKKIFERDWWSSHSKSLKILE
MW	~48kDa
Uniprot ID	Q5IXY0
Solubility (25°C)	2M Urea, PH7.4
Expiration Date	6 months from date of receipt.