

Product Datasheet

D14L Recombinant Protein (orb1495099)

Description	D14L Recombinant Protein
Species/Host	E. coli
Preservatives	2M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag;Sumo-Tag
Note	For research use only
Protein Sequence	HMKVESVTFLTLGIGCVLSYCTIPSRPINMKFKNSVETDANYNIGDTIEYLCLPGY RKQKMGPIYAKCTGTGWTLFNQCIKRRCPSPRDIDNGQLDIGGVDFGSSITYSCN SGYHLIGESKSYCELGSTGSMVWNPEAPICESVKCQSPPSISNGRHNGYEDFYID GSIVTYSCNSGYSLIGNSGVMCSGGEWSNPPTCQIVKCPHPISNGKLLAALE
MW	~37kDa
Uniprot ID	Q77HN7
Solubility (25°C)	2M Urea, PH7.4
Expiration Date	6 months from date of receipt.