

Product Datasheet

D13L Recombinant Protein (orb1495101)

Description	D13L Recombinant Protein
Species/Host	E. coli
Preservatives	0.5M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only
Protein Sequence	HMDTIKIFNHGEFDTIRNELVNLLKVVKWNTINSNVTVSSTDTIDISDCIREILYKQ FKNVRNIEVSSDISFIKYNRFNDMTITDDNMGYYLVIYLNRTKSLKTLIYPTLETVIS SSEDIMFSKSLNFRFENTKREYKLVMCISLIYKPSICRIQYDNNKYIDISDSQEGN NLCYCVITMDPHHLIDLETMCVLVDKSGKCLLVNEFYIRFRKNHIYNSFVDICMDH IFELPNTKELFTLYNDDGRNIAWDNDKLESNNTWIPKTDDEYKFLSKLMNIAKFN NTKFDYYVLVGDTDPCTVFTFKVTKYYININYELE
MW	~33kDa
Uniprot ID	Q5QC82
Solubility (25°C)	0.5M Urea, PH7.4
Expiration Date	6 months from date of receipt.