

## Product Datasheet

### TCL1A Recombinant Protein (orb1495102)

<b>Description</b>	TCL1A Recombinant Protein
<b>Species/Host</b>	E. coli
<b>Preservatives</b>	0.5M Urea, pH 7.4
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Tag</b>	His-tag
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MAECPTLGEAVTDHPDRLWAWKEFVYLDEKQHAWLPLTIEIKDRLQLRVLLRRE DVVLGRPMTPTQIGPSLLPIMWQLYPDGRYRSDSSFWRLVYHIKIDGVEDMLLE LLPDD
<b>MW</b>	~17kDa
<b>Uniprot ID</b>	<b>P56279</b>
<b>Solubility (25°C)</b>	0.5M Urea, PH7.4
<b>Expiration Date</b>	6 months from date of receipt.