

## Product Datasheet

### A33R Recombinant Protein (orb1495105)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | A33R Recombinant Protein   |
| <b>Species/Host</b>      | E. coli  |
| <b>Preservatives</b>     | PBS, 4M Urea, pH 7.4   |
| <b>Storage</b>           | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.   |
| <b>Tag</b>               | His-tag  |
| <b>Note</b>              | For research use only  |
| <b>Protein Sequence</b>  | HMASILNTRLRFLEKTSFYNCNDSITKEKIKIKHKGMSFVFYKPKHSTVVKYLSGGGI<br>YHDDLVLGKVTINDLKMMLFYMDLSYHGVTSSGAIYKLGSSIDRLSLNRTIVTKV<br>NNNYNNYNNYNNYNNYNNYNNYNNYNNYNNYNNYDDTFFDDDDLE |
| <b>MW</b>                | ~16kDa   |
| <b>Uniprot ID</b>        | <a href="#">Q318L8</a>   |
| <b>Solubility (25°C)</b> | PBS, 4M Urea, PH7.4  |
| <b>Expiration Date</b>   | 6 months from date of receipt.   |