

Product Datasheet

A14L Recombinant Protein (orb1495107)

Description	A14L Recombinant Protein
Species/Host	E. coli
Preservatives	PBS, 4M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag;Sumo-Tag
Note	For research use only
Protein Sequence	HMNKIKNPQNPNSPNLNSPPPETRNTKVFVNNLEKDHISLYNLVKSSALE
MW	~19kDa
Uniprot ID	Q910E9
Solubility (25°C)	PBS, 4M Urea, PH7.4
Expiration Date	6 months from date of receipt.