

Product Datasheet

CD297 Recombinant Protein (orb1495132)

| | |
|--------------------------|--|
| Description | CD297 Recombinant Protein |
| Species/Host | E. coli |
| Preservatives | PBS, 4M Urea, pH 7.4 |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Tag | His-tag |
| Note | For research use only |
| Protein Sequence | SEVAIKIDFD FAPGSFDDQYQGCSKQVMEKLTQGDYFTKDIEAQKNYFRMWQKA HLAWLNQGGKVL PQNM TTTTHAVAILFYTLNSNVHSD FTRAMASVARTPQQYERSF HFKYLHYL TSAIQLLRKDSIMENGLCYEVHYRTKDVHFNAYTGATIRFGQFLST SLLKEEAQEFGNQTLFTIFTCLGAPVQYFSLKKEVLIPPYELFKVINMSYHPRGNWL QLRSTGNLSTYNCQLLKA |
| MW | ~30kDa |
| Uniprot ID | Q93070 |
| Solubility (25°C) | PBS, 4M Urea, PH7.4 |
| Expiration Date | 6 months from date of receipt. |